

EHRLICH REAGENT INSTRUCTIONS

1. Carefully open the WIM Scientific Laboratories Ehrlich Reagent's factory seal.
2. Using the provided mini tester spoon, place at least .010 to .005 Grams (**Tiny amount**) of the questionable substance into the empty testing vial. *If you are testing blotter paper, make sure you are using at least half of a tab. If you are testing for psilocybin, make sure the sample is ground up and dry. Only use the vial for powdered substances such as DMT, use a white ceramic surface for testing blotter or psilocybin.*
3. Add one or two drops of the Ehrlich Reagent into the testing vial.
4. Refer to the **color chart** (on back) to determine what is present in the sample.
5. After successfully testing your substance, mini testing spoon and testing vial will need to be completely cleaned and dried before your next use.
6. After testing, the Ehrlich Reagent bottle cap should be closed tightly and placed back into the bag to ensure no leakage or unwanted exposure occurs.
7. Also included are glow sticks and wristbands...because we love you. They may come in handy! Keep out of direct sunlight and hot temperatures (Above 120 degrees) for best results and lasting usage. Ehrlich Reagent will be good for at least 3-6 months depending on proper storage.
If reagent comes into contact with skin or eyes, wash promptly with soap and water. Latex gloves are highly recommended when handling kits.
If irritation persists, receive medical attention immediately. KEEP OUT OF REACH OF CHILDREN.



Only The Highest Quality Lab Grade Chemicals Used:
ISO 9001:2008 Quality Lab Material Certified



ENRICH COLOR REACTIONS

(Read all instructions on back first)*

SUBSTANCES:

Psilocin

Psilocybin

4-ACO-DMT

LSD

DMT

AMT

5-MEO-AMT

5-MEO-MIPT

5-MEO-DMT


0 min.

1

5

10

15



* It may take up to 15 minutes to see a color change depending on the potency of the ndole. For best results, use a white ceramic surface for testing. NBOMes and other research chemicals will have no reaction or color change.

WIM Scientific Laboratories LLC™

info@wimscilabs.com

